



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in wooden case, polished stainless steel, F1: 1 g - 5 kg, F1

€2,040.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: single weights brass nickel-plated or stainless steel polished, milligram weights
- Material case: wood padded. Milligram weights (1 mg - 500 mg) in removable plastic case

Additional Information

No.	KP-3883
Manufacturer (Brand)	KERN
EAN	4058072135706
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	5000 g
Weight min	1 g
Manufacturer	CORE

