



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN set of weights in plastic case brass fine turned, 1 g - 500 g, M1

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- Material test weight: single weights brass finely turned or stainless steel finely turned, milligram weights aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

Additional Information

No.	KP-2323
Manufacturer (Brand)	KERN
EAN	4058072130206
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Weight max	500 g
Weight min	1 g
Manufacturer	CORE

