



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN set of weights in plastic case brass fine turned, 1 g - 50 g, M1

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Material test weight: single weights brass finely turned or stainless steel finely turned, milligram weights aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

## Additional Information

---

No.	KP-2582
Manufacturer (Brand)	KERN
EAN	4058072131111
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Weight max	50 g
Weight min	1 g
Manufacturer	CORE

