



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in plastic case stainless steel fine turned, 1 mg - 100000 mg, M1

€305.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: single weights brass finely turned or stainless steel finely turned, milligram weights aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

Additional Information

No.	KP-2604
Manufacturer (Brand)	KERN
EAN	4058072131203
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	100000 mg
Weight min	1 mg
Manufacturer	CORE

