



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in plastic case stainless steel fine turned, 1 g - 2000 g, F2

**€890.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: single weights stainless steel fine turned, milligram weights aluminum (1 mg - 5 mg) / nickel silver (10 mg - 500 mg)
- Material case: plastic padded. Milligram weights (1 mg - 500 mg) in removable plastic case

## Additional Information

---

No.	KP-2757
Manufacturer (Brand)	KERN
EAN	4058072131944
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	2000 g
Weight min	1 g
Manufacturer	CORE

