



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in plastic case stainless steel fine turned, 1 g - 500 g, F2

**€470.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: single weights stainless steel fine turned, milligram weights aluminum (1 mg - 5 mg) / nickel silver (10 mg - 500 mg)
- Material case: plastic padded. Milligram weights (1 mg - 500 mg) in removable plastic case

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | KP-2914          |
| Manufacturer (Brand)  | KERN             |
| EAN                   | 4058072132552    |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes              |
| Color                 | silver           |
| Material              | Stainless steel  |
| sterile               | No               |
| Weight max            | 500 g            |
| Weight min            | 1 g              |
| Manufacturer          | CORE             |

