



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set plastic case, polished stainless steel, F1 1 mg - 2 kg, F1

€1,640.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: individual weights nickel-plated brass or polished stainless steel, milligram weights
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

Additional Information

No.	KP-1241
Manufacturer (Brand)	KERN
EAN	4058072126896
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	2000000 mg
Weight min	1 mg
Manufacturer	CORE

