



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set,, stainless steel polished E2 1 g - 500 g weight set, E2

**€750.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: polished stainless steel
- Material case: plastic

## Additional Information

---

No.	KP-1691
Manufacturer (Brand)	KERN
EAN	4058072128470
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	500 g
Weight min	1 g
Manufacturer	CORE

