



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set,, stainless steel polished E2 1 g - 200 g weight set, E2

€650.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: polished stainless steel
- Material case: plastic

Additional Information

No.	KP-1959
Manufacturer (Brand)	KERN
EAN	4058072129279
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	200 g
Weight min	1 g
Manufacturer	CORE

