



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set plastic case, polished stainless steel, F1 1 mg - 200 g, F1

**€820.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: individual weights nickel-plated brass or polished stainless steel, milligram weights
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

## Additional Information

---

No.	KP-3095
Manufacturer (Brand)	KERN
EAN	4058072133009
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	200000 mg
Weight min	1 mg
Manufacturer	CORE

