



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in plastic case stainless steel polished, 1 g - 2000 g, E2

€2,250.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: single weights stainless steel polished, milligram weights
- Material case: plastic. Milligram weights (1 mg - 500 mg) in removable plastic case
- The table shows the smallest and the largest weight in each case

Additional Information

No.	KP-3597
Manufacturer (Brand)	KERN
EAN	4058072134945
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	2000 g
Weight min	1 g
Manufacturer	CORE

