



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set plastic case, polished stainless steel, F1 1 mg - 1 kg, F1

€1,240.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: individual weights nickel-plated brass or polished stainless steel, milligram weights
- Material case: plastic padded. Milligram weights (1 mg-500 mg) in removable plastic case

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-3840 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072135515 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | silver |
| Material | Stainless steel |
| sterile | No |
| Weight max | 1000000 mg |
| Weight min | 1 mg |
| Manufacturer | CORE |

