



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN test weight block cast iron lacquered, 5 kg, M2

**€72.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: cast iron, lacquered

## Additional Information

---

No.	KP-0024
Manufacturer (Brand)	KERN
EAN	4058072124380
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Cast iron
sterile	No
Form	Block
Weight	5 kg
Manufacturer	CORE

