



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight wire,, polished stainless steel 1 mg, E1

€101.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: stainless steel

Additional Information

No.	KP-2889
Manufacturer (Brand)	KERN
EAN	4058072132446
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	rod-shaped. flexible
Weight	1 mg
Manufacturer	CORE

