



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight hook weight brass fine turned, 5 kg, M1

**Price on
request
plus VAT &
Shipping**



Description

- Material test weight: finely turned brass
- Material container (optionally available): Plastic padded

Additional Information

No.	KP-1747
Manufacturer (Brand)	KERN
EAN	4058072128593
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Form	Hook
Weight	5 kg
Manufacturer	CORE

