



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob brass fine turned, 1 kg, M3

**Price on
request
plus VAT &
Shipping**



Description

- Material test weight: finely turned brass
- Material container (optionally available): Plastic padded

Additional Information

No.	KP-0175
Manufacturer (Brand)	KERN
VGKL number	645836751
EAN	4058072124809
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Form	Button
Weight	1 kg
Manufacturer	CORE

