



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob stainless steel polished, 20 kg, E2

€2,330.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: polished stainless steel
- **Note:** Shipping weights without a protective case can very likely result in damage to the weights that exceeds the respective permissible error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases that prevent damage and keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

Additional Information

No.	KP-1434
Manufacturer (Brand)	KERN
VGKL number	645831715
EAN	4058072127732
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	Button
Weight	20 kg
Manufacturer	CORE

