



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob stainless steel polished, 5 g, E2

€63.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: polished stainless steel
- **Note:** Shipping weights without a protective case can very likely result in damage to the weights that exceeds the respective permissible error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases to prevent damage and to keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-1467 |
| Manufacturer (Brand) | KERN |
| VGKL number | 645831703 |
| EAN | 4058072127923 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | silver |
| Material | Stainless steel |
| sterile | No |
| Form | Button |
| Weight | 5 g |
| Manufacturer | CORE |

