



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob, polished stainless steel 100 g, E1

€200.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: polished stainless steel
- Material case (optionally available): Wood padded

Additional Information

No.	KP-1759
Manufacturer (Brand)	KERN
EAN	4058072128708
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	Button
Weight	100 g
Manufacturer	CORE

