



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob cast iron lacquered, 1 kg, M3

€47.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: cast iron, lacquered

Additional Information

No.	KP-1917
Manufacturer (Brand)	KERN
EAN	4058072123185
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Cast iron
sterile	No
Form	Button
Weight	1 kg
Manufacturer	CORE

