



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob brass fine turned, 20 g, M1

**Price on
request
plus VAT &
Shipping**



Description

- Material test weight: finely turned brass or finely turned stainless steel

Additional Information

No.	KP-2177
Manufacturer (Brand)	KERN
VGKL number	645834745
EAN	4058072129804
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Form	Button
Weight	20 g
Manufacturer	CORE

