




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN test weight knob stainless steel polished, 2 g, E2

**€60.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: polished stainless steel
- **Note:** Shipping weights without a protective case can very likely result in damage to the weights that exceeds the respective permissible error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases that prevent damage and keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

## Additional Information

---

No.	KP-2259
Manufacturer (Brand)	KERN
VGKL number	645831702
EAN	4058072130152
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	Button
Weight	2 g
Manufacturer	CORE

