



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob brass fine turned, 2 g, M1

€7.50
plus VAT &
Shipping



Description

- Material test weight: finely turned brass or finely turned stainless steel

Additional Information

No.	KP-2823
Manufacturer (Brand)	KERN
VGKL number	645834742
EAN	4058072132224
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Form	Button
Weight	2 g
Manufacturer	CORE

