



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN test weight knob stainless steel fine turned, 500 g, F2

**€71.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: stainless steel fine turned

## Additional Information

---

No.	KP-2829
Manufacturer (Brand)	KERN
VGKL number	645833769
EAN	4058072132255
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	Button
Weight	500 g
Manufacturer	CORE

