




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN test weight knob stainless steel , 100 g, M1

**€19.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: finely turned brass or finely turned stainless steel

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | KP-2921          |
| Manufacturer (Brand)  | KERN             |
| VGKL number           | 645834707        |
| EAN                   | 4058072123369    |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes              |
| Color                 | silver           |
| Material              | Stainless steel  |
| sterile               | No               |
| Form                  | Button           |
| Weight                | 100 g            |
| Manufacturer          | CORE             |

