



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight knob stainless steel fine turned, 1 g, M1

€8.40
plus VAT &
Shipping

Product Images



Description

- Material test weight: finely turned brass or finely turned stainless steel

Additional Information

No.	KP-2975
Manufacturer (Brand)	KERN
EAN	4058072132880
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Form	Button
Weight	1 g
Manufacturer	CORE

