



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN test weight knob brass fine turned, 100 g, M3

**€14.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: finely turned brass
- Material container (optionally available): Plastic padded

## Additional Information

---

No.	KP-3446
Manufacturer (Brand)	KERN
EAN	4058072134358
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Form	Button
Weight	100 g
Manufacturer	CORE

