



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight platelet E2: 500 mg nickel silver, E2

€29.00
plus VAT &
Shipping

Product Images



Description

- Test weight material: aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- **Note:** Shipping weights without a protective case can most likely result in damage to the weights that exceeds the respective allowable error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases to prevent damage and to keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

Additional Information

No.	KP-0772
Manufacturer (Brand)	KERN
VGKL number	645831809
EAN	4058072126612
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Nickel silver
sterile	No
Form	Platelet
Weight	500 mg
Manufacturer	CORE

