



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN test weight platelet aluminum, 5 mg, E2

€29.00
plus VAT &
Shipping

Product Images



Description

- Test weight material: aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- **Note:** Shipping weights without a protective case can most likely result in damage to the weights that exceeds the respective allowable error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases that prevent damage and keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

Additional Information

No.	KP-1410
Manufacturer (Brand)	KERN
VGKL number	645831803
EAN	4058072127596
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Aluminum
sterile	No
Form	Platelet
Weight	5 mg
Manufacturer	CORE

