



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-24-11-AT, aerosol-tight, incl. rotor lid, stainless steel

€1,266.10
plus VAT &
Shipping

Product Images



Description

Rotor S-24-11-AT, aerosol-tight, incl. rotor lid, stainless steel for: Centrifuge 5430/5430 R, Centrifuge 5430/5430 R - Max. Speed: 16,049 × g (12,700 rpm) - Swing-out rotor for 24 × 1.5/2.0 mL tubes - Aerosol-tight with Eppendorf QuickLock® lid - Perfect for phase separation applications such as phenol-chloroform extraction - Compatible with Centrifuges 5430 and 5430R with serial numbers from 16,000. Software update required for other serial numbers. Please contact your responsible sales partner.

Additional Information

No.	VB-1946
Manufacturer (Brand)	Eppendorf
VGKL number	106927757
EAN	4058072419905
Gross weight	2
HSNumber	84219100
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	1.80000 kg
Material	Stainless steel
Volume MAX	2000 µl
Country of origin	Germany
Packing width	0.24000 m
Packing height	0.16000 m
Packing depth	0.24000 m
Packaging volume	10000.000000 ccm
Customs tariff number	84219100
Fixed angle rotor	90°
Type rotor	Swing-out rotor
centrifugable up to	16049 g

