



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf swing-out rotor f. 8 x 15 ml tubes, for E-1901

€316.80
plus VAT &
Shipping

Product Images



Description

Swing-out rotor type A-8-17, Eppendorf

for 8 x 15 ml tubes, autoclavable.
Max. Speed: 4,400 rpm

Additional Information

No.	VB-2012
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1902
EAN	4058072420505
Transport temperature	Room temperature
autoclavable at 121°C	No
Rotor(s) is/are suitable for	Centrifuge tubes. 8 x 15 ml
Type rotor	Swing-out rotor
for model centrifuge	Eppendorf. 5702 Eppendorf. 5702 R Eppendorf. 5702 RH

