



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific BIOLiner™ swing-out rotor

€407.00
plus VAT &
Shipping

Product Images

Description

BIOLiner™ Swinging Bucket Centrifuge samples at 3,500 rpm with the Thermo Scientific™ BIOLiner™ Swinging Bucket. Description: Click Biocontainment Lid (Set of 4) Type: Swinging Bucket Capacity: 4 x 750 ml Max. RZB: 2739 x G Max. Speed: 3,500 rpm Radius (metric): 20 cm Min. radius (metric): 10.5 cm K-factor: 13307 Max. Tubing length: 124 mm For use with (devices): Thermo Scientific Sorvall Legend XT/XF, Sorvall ST 40, Multifuge X3/X3F, Megafuge 40 and SL 40/40F Centrifuges Angle: 90° Certification/Compliance: Biocontainment certified by CAMR™, HPA, Porton Down, U.K. Biocontainment: Yes Unit Size: Each

Additional Information

No.	LD-2043
Manufacturer (Brand)	Thermo Scientific
EAN	4066292016965
Transport temperature	Room temperature

