



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf DeepWell rotor incl. 2 DeepWell plate carriers

€1,868.00
plus VAT &
Shipping

Product Images



Description

DeepWell rotor type A-2-DWP, Eppendorf

incl. 2 DeepWell plate carriers, max. loading height 85 mm.
Max. 3,700 rpm

Additional Information

No.	VB-2087
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-1898
EAN	4058072421243
Transport temperature	Room temperature
autoclavable at 121°C	No
Rotor(s) is/are suitable for	Deep well plates. 2x
Type rotor	Swing-out rotor for deep well plates
for model centrifuge	Eppendorf. 5804 Eppendorf. 5804 R

