




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf swing-bucket rotor A-4-38, incl 4 round beakers à 85/100ml

**€1,119.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor A-4-38 with round buckets, 8 places, max. Ø 38 mm, incl. 4 round buckets, 100 mL for: Centrifuge 5702/5702 R/5702 RH - Maximum speed: 3,000 × g (4.400 rpm) - Standard rotor for centrifugation of blood collection tubes or conical 15-/50-mL tubes - Round buckets for tubes with 1.5 - 100 mL - Large variety of adapters for all blood collection tubes and cell culture tube types - Aerosol tightness<sup>1)</sup> Caps for safe centrifugation of hazardous samples optionally available<sup>1)</sup> Aerosol tightness tested by Public Health England, Porton Down (UK).

## Additional Information

---

No.	VB-2029
Manufacturer (Brand)	Eppendorf
EAN	4058072420673
Transport temperature	Room temperature

