




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-4-104, incl. 4 × 750 mL round beakers

€3,222.00
plus VAT &
Shipping

Product Images



Description

Rotor S-4-104, incl. 4 × 750 mL round buckets for: Centrifuge 5810/5810 R - Swing-out rotor with high capacity and a max. speed of $3,214 \times g$ (3.900 rpm)¹⁾ - Wide range of adapters offers high flexibility for tubes and bottles from 0.2 mL to 750 mL as well as plates - Aerodynamically designed round buckets for up to 36 % energy savings compared to rotor A-4-81 - Mixed loading with 2 plate buckets and 2 round buckets offers more flexibility - Aerosol-tight QuickLock caps for safe centrifugation of hazardous samples optionally available²⁾ - Buckets, Caps and adapters are autoclavable (121 °C, 20 min) ¹⁾ Please note that older centrifuges require a software upgrade when using this rotor. For further information, please contact your local Eppendorf organization or your local Eppendorf distributor. ²⁾ Aerosol tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2217
Manufacturer (Brand)	Eppendorf
VGKL number	106974000
EAN	4058072422547
Gross weight	8.8
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	8.50000 kg
Volume MAX	750 ml
Volume MIN	0.2 ml
Country of origin	Germany
Packing width	0.42000 m
Packing height	0.24000 m
Packing depth	0.42000 m
Packaging volume	43000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	3214 g

