



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor A-4-38, 8 places, max.Ø 38 mm, without buckets

€611.00
plus VAT &
Shipping

Product Images



Description

Rotor A-4-38, 8 places, max.Ø 38 mm, without buckets for: Centrifuge 5702/5702 R/5702 RH - maximum speed: 2,750 × g (4,400 rpm) - rectangular buckets for tubes with 5 - 25 mL

Additional Information

No.	VB-2040
Manufacturer (Brand)	Eppendorf
VGKL number	556940004
EAN	4058072420772
Gross weight	1.3
HSNumber	84219100
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	1.17000 kg
Volume MAX	25000 µL
Country of origin	Germany
Packing width	0.25000 m
Packing height	0.22000 m
Packing depth	0.25000 m
Packaging volume	14000.000000 ccm
Customs tariff number	84219100
centrifugable up to	3000 g

