



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor A-4-62-MTP, with 4 MTP cups

€1,942.00
plus VAT &
Shipping

Product Images



Description

Rotor A-4-62-MTP, with 4 MTP cups for: Centrifuge 5810/5810 R – Swing-out rotor with a max. speed of 3,220 × g (4,000 rpm) – Rotor A-4-62 holds tubes and bottles from 0.2 mL to 250 mL as well as plates – Aerosol-tight caps) for safe centrifugation of hazardous samples optionally available – Rotor lifetime: max. 40,000 cycles – Rotor, cups, caps and adapters are autoclavable 1) Aerosol-tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2144
Manufacturer (Brand)	Eppendorf
VGKL number	556911002
EAN	4058072421816
Gross weight	6.505
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	5.85500 kg
Country of origin	Germany
Packing width	0.31000 m
Packing height	0.22000 m
Packing depth	0.31000 m
Packaging volume	22000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	2750 g

