



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor A-4-81, for 500 mL rectangular beakers or MTP/Flex beakers

€2,977.00
plus VAT &
Shipping

Product Images



Description

Rotor A-4-81, for 500 mL rectangular buckets or MTP/Flex buckets, incl. 4 × 500 mL rectangular buckets for: Centrifuge 5810/5810 R – swing-out rotor with a max. speed of 3,220 × g (4.000 rpm) – Rotor A-4-81 holds tubes and bottles from 0.2 mL to 500 mL as well as plates – Aerosol-tight caps) for safe centrifugation of hazardous samples optionally available – Mixed loading with 2 plate buckets and 2 rectangular buckets offers more flexibility – Buckets, caps and adapters are autoclavable (121 °C, 20 min) 1) Aerosol tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2147
Manufacturer (Brand)	Eppendorf
VGKL number	556918007
EAN	4058072421847
Gross weight	8.5
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	7.65000 kg
Country of origin	Germany
Packing width	0.30000 m
Packing height	0.20000 m
Packing depth	0.30000 m
Packaging volume	18000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	3220 g

