



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor A-4-81-MTP/Flex, swing-out rotor, incl. 4 MTP/Flex beakers

€2,977.00
plus VAT &
Shipping

Product Images



Description

Rotor A-4-81-MTP/Flex, swing-out rotor, incl. 4 MTP/Flex buckets for: Centrifuge 5810/5810 R - Swing-out rotor with a max. speed of 3,220 × g (4,000 rpm) - Rotor A-4-81 holds tubes and bottles from 0.2 mL to 500 mL as well as plates - Aerosol-tight caps¹⁾ for safe centrifugation of hazardous samples optionally available - Mixed loading with 2 plate buckets and 2 round buckets offers more flexibility - Rotor lifetime: max. 100.000 cycles - Beakers, caps and adapters are autoclavable¹⁾ Aerosol tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2154
Manufacturer (Brand)	Eppendorf
VGKL number	556925003
EAN	4058072421915
Gross weight	8.6
Transport temperature	Room temperature
HSNumber	84219100
Storage temperature	Room temperature
Net weight	7.74000 kg
Country of origin	Germany
Packing width	0.32000 m
Packing height	0.21000 m
Packing depth	0.32000 m
Packaging volume	22000.000000 ccm
Customs tariff number	84219100
Type rotor	Swing-out rotor
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	2900 g

