



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor A-4-44, incl. 4 beakers for 2 × 50 mL conical tubes

€1,771.00
plus VAT &
Shipping

Product Images



Description

Rotor A-4-44, incl. 4 buckets for 2 × 50 mL conical tubes for: Centrifuge 5804/5804 R, Centrifuge 5810/5810 R - High-speed swing-out rotor for applications up to 4,500 × g (5,000 rpm) - Rotor A-4-44 holds tubes from 0.2 mL to 100 mL - Aerosol-tight caps¹⁾ for safe centrifugation of hazardous samples optionally available - Rotor lifetime: max. 34.000 cycles - Rotor, buckets, caps and adapters are autoclavable ¹⁾ Aerosol tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2073
Manufacturer (Brand)	Eppendorf
VGKL number	106971900
EAN	4058072421106
Gross weight	4.3
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	3.87000 kg
Country of origin	Germany
Packing width	0.32000 m
Packing height	0.23000 m
Packing depth	0.32000 m
Packaging volume	24000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	4500 g

