



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-4-72, incl. 4 × 250 mL round beakers

€2,271.00
plus VAT &
Shipping

Product Images



Description

Rotor S-4-72, incl. 4 × 250 mL round buckets for: Centrifuge 5804/5804 R, Centrifuge 5810/5810 R - Swing-out rotor with medium capacity and a max. speed of $3,234 \times g$ (4.200 rpm)¹⁾ - Numerous adapters available for all common tube formats from 0.2 mL to 250 mL - Aerodynamically designed round buckets for up to 48 % energy savings compared to rotor A-4-62 - Buckets and adapters are autoclavable (121 °C, 20 min)¹⁾ Please note that older centrifuges require a software upgrade when using this rotor. For further information, please contact your local Eppendorf organization or your local Eppendorf distributor.

Additional Information

No.	VB-2091
Manufacturer (Brand)	Eppendorf
VGKL number	106974607
EAN	4058072421281
Gross weight	5
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	4.50000 kg
Volume MAX	250 ml
Volume MIN	0.2 ml
Country of origin	Germany
Packing width	0.30000 m
Packing height	0.20000 m
Packing depth	0.30000 m
Packaging volume	18000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	3234 g

