



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-4x1000 high-capacity buckets, incl. high-capacity buckets

€4,097.45
plus VAT &
Shipping

Product Images



Description

Rotor S-4x1000 high-capacity buckets, incl. high-capacity buckets for: Centrifuge 5920 R - swing-out rotor with extra large capacity and a max. speed of $3,153 \times g$ (3,700 rpm) - buckets hold both plates and tubes - buckets and adapters are autoclavable

Additional Information

No.	VB-2278
Manufacturer (Brand)	Eppendorf
EAN	4058072423117
Gross weight	0.97
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.87500 kg
Volume MAX	250 ml
Volume MIN	0.2 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	3153 g

