



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor S-4x1000 plate/tube buckets, aerosol-tight, incl. plate/tube buckets

**€3,405.15**  
plus VAT &  
Shipping

### Product Images



## Description

Rotor S-4x1000 plate/tube buckets, aerosol-tight, incl. plate/tube buckets for: Centrifuge 5920 R - High-capacity swing-out rotor with a max. speed of  $3,076 \times g$  (3,700 rpm) - Bucket holds both plates and tubes - Aerosol-tight QuickLock caps for safe centrifugation of hazardous samples optionally available<sup>1)</sup> - Buckets and adapters are autoclavable<sup>1)</sup> Aerosol-tightness tested by Public Health England, Porton Down (UK).

## Additional Information

No.	VB-2277
Manufacturer (Brand)	Eppendorf
EAN	4058072423100
Gross weight	13
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	11.70000 kg
Volume MAX	50 ml
Volume MIN	0.2 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	3076 g

