




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-4x750, aerosol-tight, without beaker

€1,473.00
plus VAT &
Shipping

Product Images



Description

Rotor S-4x750, aerosol-tight, without buckets - Swing-out rotor with medium capacity and a max. speed of 4,800 × g (4.700 rpm) 1) - Wide range of adapters offers high flexibility for tubes and bottles from 0.2 mL to 750 mL as well as plates - Mixed loading with 2 plate buckets and 2 round buckets offers more flexibility - Aerosol-tight QuickLock caps for safe centrifugation of hazardous samples optionally available 2) - Buckets and adapters are autoclavable 1) For 100 V instrument, max. speed 4,031 × g / 4,300 rpm 2) Aerosol-tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2280
Manufacturer (Brand)	Eppendorf
EAN	4058072423131
Gross weight	5.8
HSNumber	84219100
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.00000 kg
Volume MAX	750 ml
Volume MIN	0.2 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5810/5810 R
centrifugable up to	4800 g

