



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor S-4xUniversal, incl. universal beaker

€3,818.00
plus VAT &
Shipping

Product Images



Description

Rotor S-4xUniversal, incl. universal buckets for: Centrifuge 5910 R, Centrifuge 5910 Ri - Large capacity swing-out rotor with a max. speed of $4,347 \times g$ (4,500 rpm) - Wide range of adapters offers high flexibility for tubes and bottles from 4 mL to 1.000 mL as well as plates - Universal cup design accommodates both plates and tubes - Universal adapters accommodate selected plate formats and conical tubes - Aerosol-tight QuickLock caps optionally available for safe centrifugation of hazardous samples¹⁾²⁾³⁾ - Cups, caps and adapters are autoclavable¹⁾ Aerosol tightness tested by Public Health England, Porton Down (UK). ²⁾Aerosol tight caps compatible with MTP, DWP 96 500 µL, DWP 384 200 µL. Aerosol tight cap not compatible with DWP 96 1,000 µL and 2,000 µL. ³⁾When used with aerosol tight caps, max. capacity is reduced to 20 x 50 mL (5 x 50 mL per adapter).

Additional Information

No.	VB-2299
Manufacturer (Brand)	Eppendorf
EAN	4058072423322
Gross weight	11.6
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	10.44000 kg
Volume MAX	750 ml
Volume MIN	4 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5910 R
centrifugable up to	4347 g

