



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Fixed Angle Rotor FA-45-24-11 -Kit, incl rotor lid, aerosol-tight

€889.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-45-34-11, aerosol-tight, incl. rotor lid, for 5430/5430 R

for: Centrifuge 5430/5430 R, Centrifuge 5430/5430 R

- Max. Speed: $20.817 \times g$ (14.000 rpm)
- Max. Capacity: $30 \times 1.5/2.0$ mL tubes
- Aerosol-tight for safe centrifugation of hazardous samples

Additional Information

No.	VB-1941
Manufacturer (Brand)	Eppendorf
EAN	4058072419868
Transport temperature	Room temperature

