



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf fixed-angle rotor FA-45-30-11, incl rotor lid, aerosol-tight

€692.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-45-30-11, aerosol-tight, incl. rotor lid, for 5430/5430 R
for: Centrifuge 5430/5430 R, Centrifuge 5430/5430 R

- Max. Speed: 20.817 × g (14.000 rpm)
- Max. Capacity: 30 × 1.5/2.0 mL tubes
- Aerosol-tight for safe centrifugation of hazardous samples

Additional Information

No.	VB-1942
Manufacturer (Brand)	Eppendorf
EAN	4058072170806
Transport temperature	Room temperature

