



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Fixed Angle Rotor f. 6 x 8s, 6 x 5s PCR strips or 48 x 0.2 ml

€682.00
plus VAT &
Shipping

Product Images



Description

Fixed-angle rotors, Eppendorf

autoclavable

Order No. E-1899

Fixed-angle rotor type F-45-48-PCR

for 6 x 8 PCR strips, 6 x 5 PCR strips or 48 x 0.2 ml PCR tubes, max. 12,000 rpm

Additional Information

No.	VB-2082
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1899
EAN	4058072421199
Transport temperature	Room temperature
autoclavable at 121°C	No
Rotor(s) is/are suitable for	PCR Strips (8s): 6 x 0.2 ml. PCR Strips (5s): 5 x 0.2 ml or PCR Tubes: 48 x 0.2 ml
Type rotor	Fixed angle rotor
for model centrifuge	Eppendorf. 5804 Eppendorf. 5804 R

