



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-45-12-17, aerosol-tight, incl. rotor lid, aluminum

€927.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-45-12-17, aerosol-tight, incl. rotor lid, aluminum for: Centrifuge 5427R, Centrifuge 5427 R - Max. Speed: 20.598 × g (14.000 rpm) - Max. Capacity: 12 × Eppendorf Tubes® 5.0 mL - Aerosol-tight with Eppendorf QuickLock® lid

Additional Information

No.	VB-1843
Manufacturer (Brand)	Eppendorf
EAN	4058072418922
Gross weight	2.34
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.20000 kg
Material	Aluminum
Volume MAX	5000 µl
Country of origin	Germany
Packing width	0.24000 m
Packing height	0.11000 m
Packing depth	0.24000 m
Packaging volume	7000.000000 ccm
Customs tariff number	84219100
Fixed angle rotor	45°
Type rotor	Fixed angle rotor
centrifugable up to	20598 g

