



Eppendorf Rotor F-45-32-5-PCR

Product Images



Description

Rotor F-45-32-5-PCR for: Centrifuge 5424/5424 R (DP/obsolete), Centrifuge 5424/5424 R (DP) - Max. Speed: 18,615 × g (15,000 rpm) - Max. Capacity: 4 × PCR strips

Additional Information

No.	VB-1879
Manufacturer (Brand)	Eppendorf
VGKL number	106904000
EAN	4058072419257
Gross weight	1
HSNumber	84219100
Storage temperature	Raumtemperatur
Net weight	0.90000 kg
Volume MAX	200 µl
Country of origin	Deutschland
Packing width	0.24000 m
Packing height	0.11000 m
Packing depth	0.24000 m
Packaging volume	7000.00000000 ccm
Customs tariff number	84219100
centrifugable up to	17829 g

